

MGEA5 Polyclonal Antibody

Catalog No.	A9033	Category	Polyclonal Antibodies
Applications	WB	Observed MW	120kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	76kDa/95kDa/96kDa/102kDa

Immunogen Information

Immunogen Recombinant fusion protein containing a sequence corresponding to amino acids 380-580 of human MGEA5 (NP_001135906.1).

Gene ID 10724

Swiss prot O60502

Synonyms MGEA5; MEA5; NCOAT; OGA; protein O-GlcNAcase

Product information

Source Rabbit

Isotype IgG


Purification method Affinity purification

Storage Store at -20°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:200 -
1:2000



Western blot - MGEA5 Polyclonal Antibody (A9033)