

DAAM2 Polyclonal Antibody

Catalog No.	A7463	Category	Polyclonal Antibodies
Applications	WB, IHC	Observed MW	123kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	123kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-86 of human DAAM2 (NP_056160.2).
Gene ID	23500
Swiss prot	Q86T65
Synonyms	DAAM2; dj90A20A.1; disheveled-associated activator of morphogenesis 2

Product information


Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background


Recommended Dilutions

WB 1:500 -
1:2000

IHC 1:50 -
1:200



Western blot - DAAM2 Polyclonal Antibody (A7463)



Immunohistochemistry - DAAM2 Polyclonal Antibody (A7463)