


## Anti-CD54 Antibody [ABT-CD54] (A100027)

### Specifications:

|               |  |
|---------------|--|
| Name:         | Anti-CD54 Antibody [ABT-CD54]  |
| Description:  | Mouse monoclonal (ABT-CD54) antibody to CD54.  |
| Applications: | IHC  |
| Reactivity:   | Human  |
| Immunogen:    | Synthetic peptide derived from human CD54.   |
| Host:         | Mouse  |
| Clonality:    | Monoclonal   |
| Clone ID:     | ABT-CD54   |
| Isotype:      | IgG  |
| Conjugate:    | Unconjugated   |
| Purification: | Immunogen affinity purified.   |
| Product Form: | Liquid   |
| Formulation:  | Supplied in Tris Buffered Saline with 0.5% Bovine Serum Albumin and 0.1% Kathon CG/ICP.      |
| Disclaimer:   | This product is for research use only. It is not intended for diagnostic or therapeutic use. |

### Images:



Immunohistochemical analysis of paraffin-embedded human Tonsil using Anti-CD54 Antibody.