

PRODUCT DATA SHEET

BOVINE WHOLE BLOOD, NON-STERILE, SODIUM HEPARIN

Product Codes

37132-1

Pack size:

500 mL

Description:

Healthy animals aged > 1yr, mixed gender and breed; collected on wet ice following USDA guidelines, using Sodium Heparin solution at 10 units/ml.

Physical State:

Fresh, non-frozen liquid

Testing:

No analytical testing is performed on this product. The blood is visually inspected at the time of bottling to ensure the absence of clots.

Packaging, shipping/storage:**Packaging**

Polypropylene bottle

Storage Temperature

2 – 8 °C

Shipping Conditions

Shipped on the day of collection on wet ice.

Expiration:

This product is perishable. The exact shelf life will depend upon the user's application.

P/N: 74062 Rev 01