

SERPINA12 Polyclonal Antibody

Catalog No.	A8074	Category	Polyclonal Antibodies
Applications	WB	Observed MW	47kDa
Cross-reactivity	Human, Mouse	Calculated MW	47kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-320 of human SERPINA12 (NP_776249.1).
Gene ID	145264
Swiss prot	Q8IW75
Synonyms	SERPINA12; OL-64; serpin A12


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:2000



Western blot - SERPINA12 Polyclonal Antibody (A8074)