

PKNOX1 Polyclonal Antibody

Catalog No.	A15301	Category	Polyclonal Antibodies
Applications	WB	Observed MW	48kDa
Cross-reactivity	Human, Mouse	Calculated MW	44kDa/47kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 322-436 of human PKNOX1 (NP_004562.2).
Gene ID	5316
Swiss prot	P55347
Synonyms	PKNOX1; PREP1; pknox1c; homeobox protein PKNOX1


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:200 -
1:2000



Western blot - PKNOX1 Polyclonal Antibody (A15301)