

FKBP10 Polyclonal Antibody

Catalog No.	A12931	Category	Polyclonal Antibodies
Applications	WB, IHC	Observed MW	70kDa
Cross-reactivity	Human, Mouse	Calculated MW	38kDa/64kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 303-582 of human FKBP10 (NP_068758.3).
Gene ID	60681
Swiss prot	Q96AY3
Synonyms	FKBP10; BRKS1; FKBP65; OI11; OI6; PPIASE; hFKB P65; FK506 binding protein 10

Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Background

The protein encoded by this gene belongs to the FKBP-type peptidyl-prolyl cis/trans isomerase (PPIase) family. This protein localizes to the endoplasmic reticulum and acts as a molecular chaperone. Alternatively spliced variants encoding different isoforms have been reported, but their biological validity has not been determined.

Recommended Dilutions

WB 1:500 -
1:2000

IHC 1:50 -
1:200



Western blot - FKBP10 Polyclonal Antibody (A12931)