Phone: (479) 636-4361 Fax: (479) 636-3562

Website: www.pelfreez-bio.com

Certificate of Analysis

PRODUCT: Rabbit Gamma Globulin

CODE NO: 21005

LOT NO: 30241*

DATE OF

MANUFACTURE: July 22, 2019 ORIGIN: France

STORAGE: 2 − 8 °C

HEALTHY

ANIMAL This product is made from raw materials derived from healthy animals that were

STATEMENT: found free from obvious signs of infectious disease.

PROCEDURE: Rabbit gamma globulin is purified from normal rabbit serum (either sex, healthy

animals) using different fractionation procedures established in our laboratory.

ASSAYS:

Required Assays	Specification	Result
Serum Protein	≥ 95%	>99%
Electrophoresis		
Protein by Weight, O.D.	≥ 85%	86%
280 nm E ^{1%} 1.39		
Moisture Content	For information only	3%
Bioburden	≤ 1000 cfu/g	145 cfu/g

QUALITY ASSURANCE

Quality Control Release Date: November 12, 2019

Doug Stocker

Quality Systems Technician

Meets or exceeds specifications ☑ Does not meet specifications ☐

SUITABLE FOR FURTHER MANUFACTURE OR RESEARCH PURPOSES