

PSG1 Polyclonal Antibody

Catalog No.	A6399	Category	Polyclonal Antibodies
Applications	WB, IHC	Observed MW	75kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	46kDa/47kDa/48kDa


Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 180-419 of human PSG1 (NP_001171754.1).
Gene ID	5669
Swiss prot	P11464
Synonyms	PSG1; B1G1; CD66f; DHFRP2; FL-NCA-1/2; PBG1; P S-beta-C/D; PS-beta-G-1; PSBG-1; PSBG1; PSG95; PSGGA; PSGIIA; SP1; pregnancy-specific beta-1 -glycoprotein 1


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background



Immunohistochemistry - PSG1 Polyclonal Antibody (A6399)




Immunofluorescence - PSG1 Polyclonal Antibody (A6399)


Recommended Dilutions

WB 1:500 -
1:2000

IHC 1:50 -
1:100



Western blot - PSG1 Polyclonal Antibody (A6399)



Immunohistochemistry - PSG1 Polyclonal Antibody (A6399)