

[KO Validated] ARID2 Polyclonal Antibody

Catalog No.	A8601	Category	Polyclonal Antibodies
Applications	WB	Observed MW	250kDa
Cross-reactivity	Human	Calculated MW	191kDa/197kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1586-1835 of human ARID2 (NP_689854.2).
Gene ID	196528
Swiss prot	Q68CP9
Synonyms	ARID2; BAF200; p200; AT-rich interaction domain 2


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Background

Recommended Dilutions

WB 1:500 -
1:2000



Western blot - [KO Validated] ARID2 Polyclonal Antibody (A8601)



Western blot - [KO Validated] ARID2 Polyclonal Antibody (A8601)