


Anti-p60 Antibody [E17-L] (A8244)

Specifications:

Name:	Anti-p60 Antibody [E17-L]
Description:	Rabbit monoclonal (E17-L) antibody to p60.
Specificity:	This antibody recognizes the epitope between Ile527 - Leu542.
Applications:	IHC-P, IHC-Fr
Recommended Dilutions:	IHC: 1:100-1:200.
Reactivity:	Human
Immunogen:	Peptide derived from the C-terminal region of human p60 sequence.
Host:	Rabbit
Clonality:	Monoclonal
Clone ID:	E17-L
Isotype:	IgG
Conjugate:	Unconjugated
Product Form:	Liquid
Formulation:	Supplied in Tris Buffered Saline, pH 8.0, with 20 mg/ml BSA and 0.05% Sodium Azide.
Storage:	Store at 4°C.
Disclaimer:	This product is for research use only. It is not intended for diagnostic or therapeutic use.


Images:



Nuclear expression of p60 protein in lymphocytes of the germinal centre of the lymph node, detected with Anti-p60 Antibody (A8244). Formalin fixed, paraffin embedded human tissues (4µm sections) stained.

Anti-p60 Antibody [E17-L] (A8244)

Images continued:



Nuclear expression of p60 protein in colorectal adenocarcinoma, detected with Anti-p60 Antibody (A8244). Formalin fixed, paraffin embedded human tissues (4 μ m sections) stained.