

CX3CL1 Polyclonal Antibody

| | | | |
|-------------------------|------------|----------------------|-----------------------|
| Catalog No. | A14198 | Category | Polyclonal Antibodies |
| Applications | WB | Observed MW | 100kDa |
| Cross-reactivity | Human, Rat | Calculated MW | 42kDa |

Immunogen Information

| | |
|-------------------|---|
| Immunogen | Recombinant fusion protein containing a sequence corresponding to amino acids 40-180 of human CX3CL1 (NP_002987.1). |
| Gene ID | 6376 |
| Swiss prot | P78423 |
| Synonyms | CX3CL1; ABCD-3; C3Xkine; CXC3; CXC3C; NTN; NTT ; SCYD1; fractalkine; neurotactin |


Product information

| | |
|----------------------------|--|
| Source | Rabbit |
| Isotype | IgG |
| Purification method | Affinity purification |
| Storage | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3. |

Background

Recommended Dilutions

WB 1:200 -
1:2000



Western blot - CX3CL1 Polyclonal Antibody (A14198)