

Phospho-CFL1-Y139 pAb

Catalog No.	AP0330	Category	Phosphorylated Antibodies
Applications	IHC	Observed MW	
Cross-reactivity	Human, Mouse, Rat	Calculated MW	18kDa

Immunogen Information

Immunogen	A phospho specific peptide corresponding to residues surrounding Y139 of human CFL1
Gene ID	1072
Swiss prot	P23528
Synonyms	CFL1; CFL; HEL-S-15; cofilin; cofilin-1


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

IHC 1:50 -
1:100



Immunohistochemistry - Phospho-CFL1-Y139 pAb (AP0330)