

## MTMR4 Polyclonal Antibody

<b>Catalog No.</b>	A7761	<b>Category</b>	Polyclonal Antibodies
<b>Applications</b>	WB, IF	<b>Observed MW</b>	Refer to figures
<b>Cross-reactivity</b>	Human, Mouse, Rat	<b>Calculated MW</b>	133kDa

### Immunogen Information

<b>Immunogen</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 570-820 of human MTMR4 (NP_004678.3).
<b>Gene ID</b>	9110
<b>Swiss prot</b>	Q9NYA4
<b>Synonyms</b>	MTMR4; FYVE-DSP2; ZFYVE11; myotubularin-related protein 4

### Product information


<b>Source</b>	Rabbit
<b>Isotype</b>	IgG
<b>Purification method</b>	Affinity purification
<b>Storage</b>	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

### Background

### Recommended Dilutions

WB 1:500 -  
1:2000

IF 1:50 -  
1:100



Immunofluorescence - MTMR4 Polyclonal Antibody (A7761)