

## PRODUCT DATA SHEET

### RABBIT BLOOD, NON-STERILE, WITH SODIUM CITRATE

**Product Codes:**

31130 -1

**Pack size:**

500 mL

**Description:**

Rabbit whole blood is collected from individual, healthy, young rabbits (*Oryctolagus cuniculus*) Breeds: New Zealand White and/or California White that are 8 - 12 weeks of age, mixed sex at time of slaughter, into 3.8% sodium citrate at a 9:1 ratio, bottled and stored on wet ice prior to shipment.

**Physical State:**

Fresh, non-frozen liquid

**Testing:**

No analytical testing is performed on this product. The blood is visually inspected at the time of bottling to ensure the absence of clots.

**Packaging, shipping/storage:**

**Packaging**

Polypropylene bottle with lid

**Storage Temperature**

2 – 8 °C

**Shipping Conditions**

This product is shipped fresh on wet ice at 2 – 8 °C.

**Expiration**

This product is perishable. The exact shelf life will depend on the user's application.

P/N: 74008 Rev 02