

KLHL42 Polyclonal Antibody

Catalog No.	A10373	Category	Polyclonal Antibodies
Applications	IHC	Observed MW	
Cross-reactivity	Human, Mouse	Calculated MW	56kDa


Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 206-505 of human KLHL42 (NP_065833.1).
Gene ID	57542
Swiss prot	Q9P2K6
Synonyms	KLHL42; Ctb9; KLHDC5; kelch-like protein 42


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background




Immunohistochemistry - KLHL42 Polyclonal Antibody (A10373)




Immunohistochemistry - KLHL42 Polyclonal Antibody (A10373)

Recommended Dilutions

IHC 1:50 -
1:200



Immunohistochemistry - KLHL42 Polyclonal Antibody (A10373)



Immunohistochemistry - KLHL42 Polyclonal Antibody (A10373)