


## Anti-SIT Antibody [SIT-01] (A86514)

### Specifications:

Name:	Anti-SIT Antibody [SIT-01]
Description:	Mouse monoclonal [SIT-01] antibody to SIT.
Specificity:	The antibody SIT-01 reacts with SHP2-interacting transmembrane adaptor protein (SIT) expressed exclusively in lymphoid organs. It weakly crossreacts with mouse SIT.
Applications:	Flow Cytometry, IP, WB
Reactivity:	Human
Immunogen:	Bacterially produced recombinant intracellular fragment of human SIT.
Host:	Mouse
Clonality:	Monoclonal
Clone ID:	SIT-01
Isotype:	IgG1
Conjugate:	Unconjugated
Purification:	Protein A chromatography.
Concentration:	1 mg/ml
Purity:	> 95% (by SDS-PAGE).
Product Form:	Liquid
Formulation:	Supplied in Phosphate Buffered Saline, pH 7.4, with 15 mM Sodium Azide.
Storage:	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Disclaimer:	This product is for research use only. It is not intended for diagnostic or therapeutic use.

### Images:



Western blot of human Jurkat T cell line lysate (1% laurylmaltoside); non-reduced sample, immunostained by Anti-SIT Antibody (A86514) and Goat Anti-Mouse IgG Antibody (H+L) (HRP).