

CNOT10 Polyclonal Antibody

Catalog No.	A12164	Category	Polyclonal Antibodies
Applications	WB	Observed MW	90kDa
Cross-reactivity	Human, Mouse	Calculated MW	49kDa/76kDa/79kDa/82kDa/89kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 465-744 of human CNOT10 (NP_056257.1).
Gene ID	25904
Swiss prot	Q9H9A5
Synonyms	CNOT10; CCR4-NOT transcription complex subunit 10


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:2000



Western blot - CNOT10 Polyclonal Antibody (A12164)