

SVOP Polyclonal Antibody

Catalog No.	A15859	Category	Polyclonal Antibodies
Applications	WB	Observed MW	61kDa
Cross-reactivity	Human	Calculated MW	60kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human SVOP (NP_061181.1).
Gene ID	55530
Swiss prot	Q8N4V2
Synonyms	SVOP; SV2 related protein


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:2000



Western blot - SVOP Polyclonal Antibody (A15859)