


## Anti-Desmin Antibody [E18-V] (A8235)

### Specifications:

Name:	Anti-Desmin Antibody [E18-V]
Description:	Rabbit monoclonal (E18-V) antibody to Desmin.
Specificity:	This antibody recognizes the epitope between Glu105 - Pro115.
Applications:	IHC-P, IHC-Fr
Recommended Dilutions:	IHC: 1:100-1:200.
Reactivity:	Human
Immunogen:	Peptide derived from the C-terminal sequence of human desmin.
Host:	Rabbit
Clonality:	Monoclonal
Clone ID:	E18-V
Isotype:	IgG
Conjugate:	Unconjugated
Product Form:	Liquid
Formulation:	Supplied in Tris Buffered Saline, pH 8.0, with 20 mg/ml BSA and 0.05% Sodium Azide.
Storage:	Store at 4°C.
Disclaimer:	This product is for research use only. It is not intended for diagnostic or therapeutic use.

### Images:



Desmin expression detected in the uterine leiomyoma tissue with Anti-Desmin Antibody (A8235). Formalin fixed, paraffin embedded human tissue (4 $\mu$ m section) stained.