

PNMA2 Polyclonal Antibody

Catalog No.	A4448	Category	Polyclonal Antibodies
Applications	WB	Observed MW	41kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	41kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human PNMA2 (NP_009188.1).
Gene ID	10687
Swiss prot	Q9UL42
Synonyms	PNMA2; MA2; MM2; RGAG2; paraneoplastic antigen Ma2


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:2000



Western blot - PNMA2 Polyclonal Antibody (A4448)