

RPP30 Polyclonal Antibody

Catalog No.	A12560	Category	Polyclonal Antibodies
Applications	WB	Observed MW	29kDa
Cross-reactivity	Mouse, Rat	Calculated MW	29kDa/35kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-268 of human RPP30 (NP_006404.1).
Gene ID	10556
Swiss prot	P78346
Synonyms	RPP30; TSG15; ribonuclease P/MRP subunit p30


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:2000



Western blot - RPP30 Polyclonal Antibody (A12560)