

FGD5 Polyclonal Antibody

Catalog No.	A15212	Category	Polyclonal Antibodies
Applications	WB	Observed MW	178kDa
Cross-reactivity	Mouse	Calculated MW	61kDa/159kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1000-1200 of human FGD5 (NP_689749.3).
Gene ID	152273
Swiss prot	Q6ZNL6
Synonyms	FGD5; ZFYVE23; FYVE, RhoGEF and PH domain containing 5


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:2000



Western blot - FGD5 Polyclonal Antibody (A15212)