

RHCG Polyclonal Antibody

Catalog No.	A4739	Category	Polyclonal Antibodies
Applications	WB	Observed MW	57kDa
Cross-reactivity	Human, Mouse	Calculated MW	53kDa

Immunogen Information

Immunogen	A synthetic peptide of human RHCG
Gene ID	51458
Swiss prot	Q9UBD6
Synonyms	RHCG; C15orf6; PDRC2; RHGK; SLC42A3; Rh family C glycoprotein

Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at 4°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, pH7.3.

Background

Recommended Dilutions

WB 1:200 -
1:2000



Western blot - RHCG Polyclonal Antibody (A4739)