

Panitumumab Biosimilar - Research Grade

Bulk Panitumumab Biosimilar - Research Grade

Product Benefits:

ichorbio's Panitumumab Biosimilar - Research Grade is manufactured in a cGMP compliant facility. ichorbio: the best antibodies for *in vivo* research.

Size:

ichorbio's research grade panitumumab biosimilar is available in the following sizes: 5mg, 10mg, 20mg, 50mg, 100mg ichorbio regularly manufactures bulk multi-gram amounts of our panitumumab biosimilar - please contact us for pricing.

Target:

EGFR

Clone:

n/a

Isotype:

IgG2 kappa

Other Names:

Epidermal growth factor receptor, Proto-oncogene c-ErbB-1, Receptor tyrosine-protein kinase erbB-1, ERBB, ERBB1, HER1

Host:

Humanized

Species Reactivity:

Human

Specificity:

Detects human EGFR.

Purification Method:

This monoclonal antibody was purified using multi-step affinity chromatography methods such as Protein A or G depending on the species and isotype.

Background:

ichor.bio // hello@ichor.bio

Products are for research use only. Not for use in diagnostic or therapeutic procedures.

ichorbio, ichorbio logo and all other trademarks are the property of ichorbio LTD © ichorbio 2023



Panitumumab biosimilar is a recombinant, human IgG2 kappa monoclonal antibody that binds specifically to the human epidermal growth factor receptor (EGFR). It has an approximate molecular weight of 147 kDa and is produced in genetically engineered mammalian (Chinese hamster ovary) cells.

Concentration:

1.0 - 5.0 mg/ml

Formulation:

Sterile, colorless liquid with a pH value of 5.6 to 6.0, containing sodium chloride and sodium acetate. BSA and Azide free.

Purity:

>95% by SDS-PAGE and HPLC

Endotoxin:

? 0.75 EU/mg as determined by the LAL method

Storage:

This antibody is stable for at least 4 weeks when stored at 2-8°C. For long term storage, aliquot in working volumes without diluting and store at – 20°C or -80°C. Avoid repeated freeze thaw cycles.

Applications:

Functional Assays

Application Notes:

Each investigator should determine their own optimal working dilution for specific applications.

Use: ichorbio's panitumumab biosimilar is for research use only (RUO). **Isotype Control:**

[Bulk Human IgG2 Isotype Control \(IB2\) \[ICH2255\]](#)