

Anti-N Cadherin Antibody [MNCD2] (A248060)

Specifications:


Name:	Anti-N Cadherin Antibody [MNCD2]
Description:	Rat monoclonal [MNCD2] antibody to N Cadherin.
Specificity:	This antibody recognizes a protein of ~140kDa, identified as N-Cadherin (NCAD), also known as CD325. NCAD is a member of the Cadherin superfamily, and consists of five extracellular repeats, a transmembrane domain and a cytoplasmic domain. CD325 deficient mice die at day 10 of gestation and embryos display major heart defects and malformed neural tubes and somites. Consistent with this, CD325 has been implicated in several aspects of cardiac development including the precardiac mesoderm, establishment of left-right symmetry and cardiac looping morphogenesis. Furthermore, CD325 is normally involved in inducing cell cycle arrest and its expression is frequently deregulated in cancer cells. Studies have linked N-cadherin to cancer metastasis by showing the aggressive tumor cells had preferentially turned on N-cadherin as opposed to E- or P-cadherin.
Applications:	Flow Cytometry, WB, IHC-P
Recommended Dilutions:	Flow Cytometry: 1-2 µg/million cells, WB: 1-2 µg/ml, IHC-P: 1-2 µg/ml
Reactivity:	Human, Mouse, Chicken
Immunogen:	Recombinant fragment, around amino acids 308-597, of mouse CDH2 protein. The epitope is the extracellular domain of N Cadherin.
Host:	Rat
Clonality:	Monoclonal
Clone ID:	MNCD2
Isotype:	IgG2a
Conjugate:	Unconjugated
Purification:	Protein A/G chromatography.
Concentration:	200 µg/ml
Product Form:	Liquid
Formulation:	Supplied in 10mM Phosphate Buffered Saline with 0.05% BSA and 0.05% Sodium Azide.
Storage:	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
General Notes:	This monoclonal antibody is also available in a different formulation without BSA and Sodium Azide - Anti-N Cadherin Antibody [MNCD2] - BSA and Azide free (A251243).

Anti-N Cadherin Antibody [MNCD2] (A248060)


Specifications continued:

Disclaimer: This product is for research use only. It is not intended for diagnostic or therapeutic use.

Images:



Immunohistochemical analysis of formalin-fixed, paraffin-embedded human heart using Anti-N Cadherin Antibody [MNCD2].



Immunohistochemical analysis of formalin-fixed, paraffin-embedded human heart using Anti-N Cadherin Antibody [MNCD2].