

PRODUCT DATA SHEET

RABBIT GAMMA GLOBULIN

Product Codes21005-1
21005-2
21005-5**Pack size:**1 g
10 g
500 g**Description:**

Rabbit gamma globulin is purified from normal rabbit serum. (Species: *Oryctolagus cuniculus*. Breeds: New Zealand White and/or California White) Healthy animals of either sex are utilized. The serum is processed using different fractionation procedures and lyophilization.

Physical State:

Off-white, fine, flaky powder

Testing:Serum Protein Electrophoresis: $\geq 95\%$ Protein by Weight: O.D. 280 nm $E^{1\%1.39}$: $\geq 85\%$ - Residual salts, lipids, etc., comprise the rest of the weight

Moisture Content: Report result

Bioburden: ≤ 1000 CFU/g**Packaging, shipping/storage:****Packaging**

1 g and 10 g: Polypropylene bottle; 500 g: HDPE container with lid

Storage Temperature

2 – 8 °C

Shipping Conditions

Ambient temperature

Expiration

Product quality is guaranteed to meet Pel-Freez Biologicals' specifications for 1 year from the date of receipt by the customer as long as the product is stored in accordance with the indicated storage conditions.

Application Notes: Reportedly useful as a stabilizer/blocker for EIA and RIA systems, IHC assays, immunoprecipitation, Western Blot and ELISA experiments, and protein standards.

P/N: 74022 Rev 03