

PRODUCT DATA SHEET

SHEEP GAMMA GLOBULIN

Product Codes:

29005-1
29005-2
29005-3

Pack size:

1 g
10 g
50 g

Description:

Sheep gamma globulin is purified from normal sheep serum or plasma using a combination of existing methods developed in our laboratory.

Physical State:

Off white, fine, flaky powder

Testing:

Serum Protein Electrophoresis: $\geq 95\%$

Protein by Weight: O.D. 280 nm $E^{1\%1.39}$: $\geq 85\%$ - Residual salts, lipids, etc., comprise the rest of the weight.

Moisture Content: FIO

Packaging, shipping/storage:

Packaging

Polypropylene bottle

Storage Temperature

2 – 8 °C

Shipping Conditions

Ambient temperature

Expiration

Product quality is guaranteed to meet Pel-Freez Biologicals' specifications for 1 year from the date of receipt by the customer as long as the product is stored in accordance with the indicated storage conditions.

Application Notes: Reportedly useful as a stabilizer/blocker for EIA and RIA systems, IHC assays, immunoprecipitation, Western Blot and ELISA experiments, and protein standards.

P/N: 74025 Rev 01