


# Anti-CD79b Antibody [IGB/2940R] - BSA and Azide free (A253980)

## Specifications:

Name:	Anti-CD79b Antibody [IGB/2940R] - BSA and Azide free
Description:	Recombinant rabbit monoclonal [IGB/2940R] antibody to CD79b.
Applications:	ELISA
Reactivity:	Human
Immunogen:	Recombinant fragment, around amino acids 29-159, of human CD79b protein. The exact sequence is proprietary.
Host:	Rabbit
Clonality:	Monoclonal
Clone ID:	IGB/2940R
Isotype:	IgG
Conjugate:	Unconjugated
Purification:	Protein A/G chromatography.
Concentration:	1 mg/ml
Product Form:	Liquid
Formulation:	Supplied in 10mM Phosphate Buffered Saline; without Sodium Azide and carrier free.
Storage:	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
General Notes:	This monoclonal antibody is also available in a different formulation with BSA and Sodium Azide - Anti-CD79b Antibody [IGB/2940R] (A250800).
Disclaimer:	This product is for research use only. It is not intended for diagnostic or therapeutic use.

## Images:



SDS-PAGE analysis of Anti-CD79b Antibody [IGB/2940R] under non-reduced and reduced conditions; showing intact IgG and intact heavy and light chains, respectively. SDS-PAGE analysis confirms the integrity and purity of the antibody.