


Anti-Fibronectin Antibody [Fn-3] - BSA and Azide free (A251773)

Specifications:

Name:	Anti-Fibronectin Antibody [Fn-3] - BSA and Azide free
Description:	Mouse monoclonal [Fn-3] antibody to Fibronectin.
Specificity:	This antibody reacts with human cellular fibronectin but not plasma fibronectin.
Applications:	Flow Cytometry, IF, IHC-P
Recommended Dilutions:	Flow Cytometry: 1-2 µg/million cells, IF: 1-2 µg/ml, IHC-P: 1-2 µg/ml
Reactivity:	Human
Immunogen:	FR5 cells, derived by SV40 transformation of human mammary epithelial cells.
Host:	Mouse
Clonality:	Monoclonal
Clone ID:	Fn-3
Isotype:	IgG1
Light Chains:	kappa
Conjugate:	Unconjugated
Purification:	Protein A/G chromatography.
Concentration:	1 mg/ml
Product Form:	Liquid
Formulation:	Supplied in 10mM Phosphate Buffered Saline; without Sodium Azide and carrier free.
Storage:	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
General Notes:	This monoclonal antibody is also available in a different formulation with BSA and Sodium Azide - Anti-Fibronectin Antibody [Fn-3] (A248591).
Disclaimer:	This product is for research use only. It is not intended for diagnostic or therapeutic use.

Anti-Fibronectin Antibody [Fn-3] - BSA and Azide free (A251773)

Images:



SDS-PAGE analysis of Anti-Fibronectin Antibody [Fn-3] under non-reduced and reduced conditions; showing intact IgG and intact heavy and light chains, respectively. SDS-PAGE analysis confirms the integrity and purity of the antibody.