


## Anti-Band 3 Antibody [Q1/156] - BSA and Azide free (A253163)

### Specifications:

Name:	Anti-Band 3 Antibody [Q1/156] - BSA and Azide free
Description:	Mouse monoclonal [Q1/156] antibody to Band 3.
Applications:	ELISA, IF, WB
Recommended Dilutions:	IF: 1-2 µg/ml, WB: 1-2 µg/ml
Reactivity:	Human
Immunogen:	Human fetal liver cells.
Host:	Mouse
Clonality:	Monoclonal
Clone ID:	Q1/156
Isotype:	IgG1
Light Chains:	kappa
Conjugate:	Unconjugated
Purification:	Protein A/G chromatography.
Concentration:	1 mg/ml
Product Form:	Liquid
Formulation:	Supplied in 10mM Phosphate Buffered Saline; without Sodium Azide and carrier free.
Storage:	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
General Notes:	This monoclonal antibody is also available in a different formulation with BSA and Sodium Azide - Anti-Band 3 Antibody [Q1/156] (A249983).
Disclaimer:	This product is for research use only. It is not intended for diagnostic or therapeutic use.

## Anti-Band 3 Antibody [Q1/156] - BSA and Azide free (A253163)

### Images:



SDS-PAGE analysis of Anti-Band 3 Antibody [Q1/156] under non-reduced and reduced conditions; showing intact IgG and intact heavy and light chains, respectively. SDS-PAGE analysis confirms the integrity and purity of the antibody.