

UBE2D4 Polyclonal Antibody

Catalog No.	A8332	Category	Polyclonal Antibodies
Applications	WB	Observed MW	17kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	16kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-147 of human UBE2D4 (NP_057067.1).
Gene ID	51619
Swiss prot	Q9Y2X8
Synonyms	UBE2D4; HBUCE1; ubiquitin-conjugating enzyme E2 D4


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:2000



Western blot - UBE2D4 Polyclonal Antibody (A8332)