

## NOTUM Polyclonal Antibody

<b>Catalog No.</b>	A16600	<b>Category</b>	Polyclonal Antibodies
<b>Applications</b>	WB	<b>Observed MW</b>	56kDa
<b>Cross-reactivity</b>	Human	<b>Calculated MW</b>	55kDa

### Immunogen Information

<b>Immunogen</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 20-320 of human NOTUM (NP_848588.3).
<b>Gene ID</b>	147111
<b>Swiss prot</b>	Q6P988
<b>Synonyms</b>	NOTUM; hNOTUM; palmitoleoyl-protein carboxyles terase NOTUM


### Product information

<b>Source</b>	Rabbit
<b>Isotype</b>	IgG
<b>Purification method</b>	Affinity purification
<b>Storage</b>	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

### Background

### Recommended Dilutions

WB 1:500 -  
1:2000



Western blot - NOTUM Polyclonal Antibody (A16600)