

PRAP1 Polyclonal Antibody

Catalog No.	A8770	Category	Polyclonal Antibodies
Applications	WB	Observed MW	17kDa
Cross-reactivity	Human	Calculated MW	10kDa/16kDa/17kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-151 of human PRAP1 (NP_660203.3).
Gene ID	118471
Swiss prot	Q96NZ9
Synonyms	PRAP1; PRO1195; UPA; proline-rich acidic prote in 1


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:200 -
1:500



Western blot - PRAP1 Polyclonal Antibody (A8770)