

ANP32A Polyclonal Antibody

Catalog No.	A5768	Category	Polyclonal Antibodies
Applications	WB, IHC, IF	Observed MW	28kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	28kDa


Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human ANP32A (NP_006296.1).
Gene ID	8125
Swiss prot	P39687
Synonyms	ANP32A; C15orf1; HPPCn; I1PP2A; LANP; MAPM; PHAP1; PHAPI; PP32; acidic nuclear phosphoprotein 32 family member A


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background



Immunohistochemistry - ANP32A Polyclonal Antibody (A5768)




Immunohistochemistry - ANP32A Polyclonal Antibody (A5768)

Recommended Dilutions


WB 1:500 -
1:2000

IHC 1:50 -
1:200


IF 1:50 -
1:200



Western blot - ANP32A Polyclonal Antibody (A5768)



Immunohistochemistry - ANP32A Polyclonal Antibody (A5768)



Immunofluorescence - ANP32A Polyclonal Antibody (A5768)