

TBRG4 Polyclonal Antibody

Catalog No.	A0276	Category	Polyclonal Antibodies
Applications	WB, IF, IP	Observed MW	70kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	58kDa/70kDa


Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 382-631 of human TBRG4 (NP_004740.2).
Gene ID	9238
Swiss prot	Q969Z0
Synonyms	TBRG4; CPR2; FASTKD4; protein TBRG4

Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background




Immunoprecipitation - TBRG4 Polyclonal Antibody (A0276)

Recommended Dilutions


WB 1:500 -
1:2000

IF 1:50 -
1:100

IP 1:50 -
1:200



Western blot - TBRG4 Polyclonal Antibody (A0276)



Immunofluorescence - TBRG4 Polyclonal Antibody (A0276)