

SLC16A4 Polyclonal Antibody

Catalog No.	A14818	Category	Polyclonal Antibodies
Applications	WB, IHC, IF	Observed MW	Refer to figures
Cross-reactivity	Human, Rat	Calculated MW	35kDa/42kDa/46kDa/49kDa/54kDa

Immunogen Information

Immunogen	A synthetic peptide corresponding to a sequence within amino acids 250-350 of human SLC16A4 (NP_004687.1).
Gene ID	9122
Swiss prot	O15374
Synonyms	SLC16A4; MCT4; MCT5; monocarboxylate transporter 5

Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Background

Recommended Dilutions


WB 1:500 -
1:2000

IHC 1:50 -
1:200


IF 1:50 -
1:200




Immunohistochemistry - SLC16A4 Polyclonal Antibody (A14818)



Immunohistochemistry - SLC16A4 Polyclonal Antibody (A14818)



Immunohistochemistry - SLC16A4 Polyclonal Antibody (A14818)



Immunofluorescence - SLC16A4 Polyclonal Antibody (A14818)