

METTL1 Polyclonal Antibody

Catalog No.	A3958	Category	Polyclonal Antibodies
Applications	WB	Observed MW	29kDa
Cross-reactivity	Human	Calculated MW	17kDa/31kDa

Immunogen Information

Immunogen	A synthetic Peptide of human METTL1
Gene ID	4234
Swiss prot	Q9UBP6
Synonyms	METTL1; C12orf1; TRM8; TRMT8; YDL201w; methyltransferase like 1

Product information


Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at 4°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, pH7.3.

Background

This gene is similar in sequence to the *S. cerevisiae* YDL201w gene. The gene product contains a conserved S-adenosylmethionine-binding motif and is inactivated by phosphorylation. Alternative splice variants encoding different protein isoforms have been described for this gene. A pseudogene has been identified on chromosome X.

Recommended Dilutions

WB 1:500 -
1:1000



Western blot - METTL1 Polyclonal Antibody (A3958)