

SHTN1 Polyclonal Antibody

Catalog No.	A15179	Category	Polyclonal Antibodies
Applications	WB	Observed MW	Refer to figures
Cross-reactivity	Mouse	Calculated MW	23kDa, 52- 73kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human SHTN1 (NP_001120683.1).
Gene ID	57698
Swiss prot	A0MZ66
Synonyms	SHTN1; KIAA1598; shootin-1

Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:2000

