

Phospho-MAPK9/MAPK10-T183 pAb

Catalog No.	AP0395	Category	Phosphorylated Antibodies
Applications	WB, IHC	Observed MW	54kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	27kDa/44kDa/48kDa/52kDa

Immunogen Information

Immunogen	A phospho specific peptide corresponding to residues surrounding T183 of human MAPK9/MAPK10
Gene ID	5601/5602
Swiss prot	P45984/P53779
Synonyms	

Product information


Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:1000

IHC 1:50 -
1:100



Western blot - Phospho-MAPK9/MAPK10-T183 pAb (AP0395)