

NPRL2 Polyclonal Antibody

Catalog No.	A5859	Category	Polyclonal Antibodies
Applications	WB, IF	Observed MW	40kDa
Cross-reactivity	Human, Mouse	Calculated MW	23kDa/43kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human NPRL2 (NP_006536.3).
Gene ID	10641
Swiss prot	Q8WTW4
Synonyms	NPRL2; FFEVF2; NPR2; NPR2L; TUSC4; NPR2-like, GATOR1 complex subunit

Product information


Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background


Recommended Dilutions

WB 1:500 -
1:2000

IF 1:50 -
1:200



Western blot - NPRL2 Polyclonal Antibody (A5859)



Immunofluorescence - NPRL2 Polyclonal Antibody (A5859)