

[KO Validated] CSK Polyclonal Antibody

Catalog No.	A0735	Category	Polyclonal Antibodies
Applications	WB, IF	Observed MW	51kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	50kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human CSK (NP_001120662.1).
Gene ID	1445
Swiss prot	P41240
Synonyms	CSK; tyrosine-protein kinase CSK

Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background


Recommended Dilutions

WB 1:500 -
1:2000

IF 1:50 -
1:200



Western blot - [KO Validated] CSK Polyclonal Antibody (A0735)



Western blot - [KO Validated] CSK Polyclonal Antibody (A0735)