

## **Tapered Tip Collection Swab (80mm Break Point)**

**Part Number: 292666**

<b>Handle Length (mm) [in]</b>	(152.4) [6.00]
<b>Handle Diameter (mm) [in]</b>	(2.4) [0.095]
<b>Flock Head Length (mm) [in]</b>	(18.42) [0.73]
<b>Breakpoint (mm) [in]</b>	(80) [3.15]
<b>Individual Sterile Package Dimensions (mm) [in]</b>	(7 x 190) [0.28 x 7.48]
<b>Handle Material</b>	Polypropylene
<b>Flock Material</b>	Polyamide (Nylon)
<b>Sterilization Method</b>	Ethylene Oxide Gas
<b>Approved Jurisdictions</b>	USA, Canada, and EU
<b>Swab Quantity per Shelf Pack</b>	100
<b>Shelf Pack Quantity per Shipper Box</b>	10
<b>Shipper Boxes per Pallet</b>	42
<b>Shipper Box Weight (kg) [lbs]</b>	(3.40) [7.5]
<b>Shipper Box Outer Dimensions (mm) [in]</b>	(505 x 400 x 241) [19-7/8 x 15-3/4 x 9-1/2]