



## Product Data Sheet

---

**Product Name:** VIP, human, porcine, rat

**Catalog Number:** AS-22872 (0.5 mg)                      Lot Number: See label on vial  
AS-22873 (1.0 mg)

**Sequence:** H-His-Ser-Asp-Ala-Val-Phe-Thr-Asp-Asn-Tyr-Thr-Arg-Leu-Arg-Lys-Gln-  
Met-Ala-Val-Lys-Lys-Tyr-Leu-Asn-Ser-Ile-Leu-Asn-NH<sub>2</sub> (3-letter code)  
HSDAVFTDNYTRLRKQMAVKKYLNSILN-NH<sub>2</sub> (1-letter code)

**Molecular Weight:** 3552.9

**Peptide Purity:** >95%

**Appearance:** Lyophilized white powder

**Peptide Reconstitution:** VIP peptide is freely soluble in water.

**Storage:** VIP peptide is shipped at ambient temperature. Upon receipt, store lyophilized peptide at -20°C or lower. Reconstituted peptide can be aliquoted and stored at -20 °C or lower.

**Reference:** Paul, S. et al. *J. Biol. Chem.* **276**, 28314 (2001).

**Related Products:**

<b>Name</b>	<b>Cat #</b>	<b>Size</b>
VIP, human, porcine, rat, FAM-labeled	AS-23614	0.5 mg
FAM-HSDAVFTDNYTRLRKQMAVKKYLNSILN-NH <sub>2</sub>		

*For Research Use Only*