

Anti-Human IgD Antibody [IA6-2] (PerCP-Cyanine 5.5) (A121968)

Specifications:

Name:	Anti-Human IgD Antibody [IA6-2] (PerCP-Cyanine 5.5)
Description:	Mouse monoclonal [IA6-2] antibody to Human IgD (PerCP-Cyanine 5.5).
Applications:	Flow Cytometry
Recommended Dilutions:	Flow Cytometry: This reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension.
Reactivity:	Human
Immunogen:	Human IgD.
Host:	Mouse
Clonality:	Monoclonal
Clone ID:	IA6-2
Isotype:	IgG2a
Light Chains:	kappa
Conjugate:	PerCP-Cyanine 5.5
Purification:	This antibody is conjugated with tandem dye PerCP-Cyanine 5.5 under optimum conditions. The conjugate is purified by size-exclusion chromatography.
Product Form:	Liquid
Formulation:	Supplied in Phosphate Buffered Saline with 15 mM Sodium Azide.
Storage:	Store in the dark at 2-8°C. Do not freeze!
Disclaimer:	This product is for research use only. It is not intended for diagnostic or therapeutic use.