

### 100mm Nasopharyngeal (NP) Specimen Collection Swab Cross Reference

Manufacturer	Part Number	Description	Microbrush Part Number
Copan	503CS01	Flexible Mini Tip Nasopharyngeal Swab, 100mm Breakpoint	292489 or 292497
Puritan	25-3317-H	Ultrafine Tip Nasopharyngeal Swab, 100mm Breakpoint	
	25-3317-U	Ultrafine Tip Nasopharyngeal Swab, 100mm Breakpoint	
	25-3320-H EMB 100MM	Revised Ultrafine Flock Swab w/Polystyrene Handle, 80mm Breakpoint	
	25-3320-U EMB 100MM	Revised Ultrafine Flock Swab w/Polystyrene Handle, 80mm Breakpoint	

### 80mm Nasopharyngeal (NP) Specimen Collection Swab Cross Reference

Manufacturer	Part Number	Description	Microbrush Part Number
Copan	525CS01	Ultra Thin Minitip Nasopharyngeal Swab, 80mm break point	292496
Puritan	25-3320-H EMB 80MM	Revised Ultrafine Flock Swab w/Polystyrene Handle, 80mm Breakpoint	
	25-3320-U EMB 80MM	Revised Ultrafine Flock Swab w/Polystyrene Handle, 80mm Breakpoint	

### Microbrush Nasopharyngeal Swab Details

PART #	COLOR	STERILE / NON-STERILE	BREAK POINT (MM)	OVERALL LENGTH (MM)	TIP LENGTH (MM)	MAX. TIP DIAMETER (MM)	PACKAGING
292489	White	Sterile	100	152	18.42	4.06	300 swabs/pack, 5 packs/case
292497	White	Sterile	100	152	18.42	4.06	100 swabs/pack, 10 packs/case
292496	White	Sterile	80	152	18.42	4.06	100 swabs/pack, 10 packs/case

For additional information or product samples, please visit <https://info.microbrush.com/test-swabs>, email [customerservice@microbrush.com](mailto:customerservice@microbrush.com), or call at 866-866-8698.