

C17orf70 Polyclonal Antibody

Catalog No.	A8591	Category	Polyclonal Antibodies
Applications	WB	Observed MW	110kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	27kDa/77kDa/93kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 652-881 of human C17orf70 (NP_079437.5).
Gene ID	80233
Swiss prot	Q0VG06
Synonyms	FAAP100; C17orf70; Fanconi anemia core complex-associated protein 100


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:2000



Western blot - C17orf70 Polyclonal Antibody (A8591)