

ENC1 Polyclonal Antibody

Catalog No.	A10257	Category	Polyclonal Antibodies
Applications	WB	Observed MW	66kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	57kDa/66kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human ENC1 (NP_003624.1).
Gene ID	8507
Swiss prot	O14682
Synonyms	ENC1; CCL28; ENC-1; KLHL35; KLHL37; NRPB; PIG1 0; TP53I10; ectodermal-neural cortex 1

Product information


Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

This gene encodes a member of the kelch-related family of actin-binding proteins. The encoded protein plays a role in the oxidative stress response as a regulator of the transcription factor Nrf2, and expression of this gene may play a role in malignant transformation. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

Recommended Dilutions

WB 1:1000 -
1:2000



Western blot - ENC1 Polyclonal Antibody (A10257)