
SDS-PAGE Running Buffer (10X) (Catalog# TBS5015)**DESCRIPTION**

10x Tris/Glycine/SDS Electrophoresis Buffer (1000 mL) is a 10x premixed protein electrophoresis buffer for SDS-PAGE electrophoresis. 1 X running buffer contains 25 mM Tris, 192 mM glycine, and 0.1% SDS, pH 8.3.

Directions for Use: Dilute Tris/Glycine/SDS Electrophoresis Buffer (10X) to a 1X running buffer using ddH₂O. This product supplies enough 10X material to make 10 liters of 1X solution.

Storage: Supplied as a 10X solution, 1000 mL. Store at room temperature.

This product is for *IN VITRO* research use only and is not intended for use in humans or animals.

This product is not intended for use as a therapeutic or in diagnostic procedures.

Please visit www.tribioscience.com for a complete listing of recommended companion products.

Orders Tel: 650 755-8332
order@tribioscience.com

Support Tel: 650 755-8332
info@tribioscience.com

Web: www.tribioscience.com

RELATED PRODUCTS:

Protein Assay Kit (catalog# TBS2005)

PBS-1x (catalog# TBS5003)

PBS-10x (catalog# TBS5004)

BSA Standard solution (catalog# TBS5002)

Cell Lysis Buffer (catalog# TBS5001)

6x Protein Loading Buffer (catalog# TBS5014)

4x PAGE-gel transfer buffer (catalog# TBS5016)

RIPA Buffer (catalog# TBS5017)