

OMD Polyclonal Antibody

Catalog No.	A12256	Category	Polyclonal Antibodies
Applications	WB	Observed MW	45kDa
Cross-reactivity	Mouse	Calculated MW	49kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 122-421 of human OMD (NP_005005.1).
Gene ID	4958
Swiss prot	Q99983
Synonyms	OMD; OSAD; SLRR2C; osteomodulin


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:2000



Western blot - OMD Polyclonal Antibody (A12256)